

```
classdef TopoDatastore < matlab.io.Datastore & ...
    matlab.io.datastore.Minibatchable & ...
    matlab.io.datastore.Shuffleable
%TOPODATASTORE Custom datastore to read in data from file structure for TMN
% To be used as a datastore to a local directory of TUEG files.
%
% Authors:
%     Michael Caiola (Michael.Caiola@fda.hhs.gov)
%     Meijun Ye (Meijun.Ye@fda.hhs.gov)
%
% Disclaimer: This software and documentation (the "Software") were
% developed at the Food and Drug Administration (FDA) by employees of
% the Federal Government in the course of their official duties.
% Pursuant to Title 17, Section 105 of the United States Code,
% this work is not subject to copyright protection and is in the
% public domain. Permission is hereby granted, free of charge, to any
% person obtaining a copy of the Software, to deal in the Software
% without restriction, including without limitation the rights to
% use, copy, modify, merge, publish, distribute, sublicense, or sell
% copies of the Software or derivatives, and to permit persons to
% whom the Software is furnished to do so. FDA assumes no
% responsibility whatsoever for use by other parties of the Software,
% its source code, documentation or compiled executables, and makes
% no guarantees, expressed or implied, about its quality,
% reliability, or any other characteristic. Further, use of this code
% in no way implies endorsement by the FDA or confers any advantage
% in regulatory decisions. Although this software can be
% redistributed and/or modified freely, we ask that any derivative
% works bear some notice that they are derived from it, and any
% modified versions bear some notice that they have been modified.
% The Software is not intended to make clinical diagnoses or to be
% used in any way to diagnose or treat subjects for whom the EEG is
% taken.

properties
    Datastore
    Labels
    NumClasses
    SequenceDimension
    MiniBatchSize
    Augmented
    TestSet
    newHz

    fband=[1 4; 4 8; 8 12; 12 16; 16 20; 25 40];

    chLocations
end

properties(SetAccess = protected)
    NumObservations
```

```
end

properties(Access = private)
    CurrentFileIndex
    %FileSet matlab.io.datastore.DsFileSet
end

methods
    function ds = TopoDatastore(folder,newHz,chLocations,varargin)
        % ds = NewDatastore(folder,newHz) creates a custom datastore
        % from the data in folder with new sampling rate newHz
        %
        % chLocations are needed for topoplot, use chlocs2.mat
        %
        % Optional arguments:
        %     "DataAugmentation" (default = false)
        %         Augments the data in random shuffles of 90 seconds
        %     "TestSet" (default = false)
        %         Creates a testset of 90 second data, to be used when
        %         DataAugmentation is being used

        if isempty(newHz)
            ds.newHz = 250;
        else
            ds.newHz = newHz;
        end

        inputs = ds.parseInputs(varargin{:});

        if inputs.DataAugmentation
            ds.Augmented = true;
        else
            ds.Augmented = false;
        end

        if inputs.TestSet
            ds.TestSet = true;
            ds.Augmented = false;
        else
            ds.TestSet = false;
        end

        if exist("topoplot",'file')~=2
            %open eeglab
            warning("Need EEGLAB! Trying to open...")
            try
                path = pwd;
                eeglabpath = uigetdir([], "Locate EEGLAB Directory");
                cd(eeglabpath)
                eeglab
                close
            catch
            end
        end
    end
end
```

```
        cd(path)
    catch
        error("Could not find EEGLAB. Open prior to use to add to path.")
    end
end

% Save chLocations
ds.chLocations = chLocations;

% Create file datastore.
fds = fileDatastore(folder, ...
    'ReadFcn',@(x) readSequence(x,ds.newHz,ds.Augmented,ds.TestSet,...
    ds.fband,ds.chLocations), ...
    'IncludeSubfolders',true);
ds.Datastore = fds;

% Read labels from folder names.
numObservations = numel(fds.Files);
for i = 1:numObservations
    file = fds.Files{i};
    filepath = fileparts(file);
    [~,label] = fileparts(filepath);
    labels{i,1} = label;
end
ds.Labels = categorical(labels);
ds.NumClasses = numel(unique(labels));

% Determine sequence dimension.
X = preview(fds);
ds.SequenceDimension = size(X,1);

% Initialize datastore properties.
ds.MinibatchSize = 128;
ds.NumObservations = numObservations;
ds.CurrentFileIndex = 1;
end

function tf = hasdata(ds)
    % tf = hasdata(ds) returns true if more data is available.

    %tf = hasdata(ds.Datastore);
    tf = ds.CurrentFileIndex + ds.MinibatchSize - 1 ...
        <= ds.NumObservations;
end

function [data,info] = read(ds)
    % [data,info] = read(ds) read one mini-batch of data.

    minibatchSize = ds.MinibatchSize;
    if ds.NumObservations < ds.MinibatchSize
        ds.MinibatchSize = ds.NumObservations;
```

```
end

i = 0;
while i < miniBatchSize && hasdata(ds) %hasdata(ds.Datastore)
    i = i + 1;
    Predictors{i,1} = read(ds.Datastore);
    %temppred = read(ds.Datastore);
    %predictors{i,1} = spect(temppred);
    Response(i,1) = ds.Labels(ds.CurrentFileIndex);
    ds.CurrentFileIndex = ds.CurrentFileIndex + 1;
end

data = preprocessData(Predictors,Response);
info.Size = size(data);
end

function reset(ds)
    % reset(ds) resets the datastore to the start of the data.

    reset(ds.Datastore);
    ds.CurrentFileIndex = 1;
end

function dsNew = subset(ds,in)
    dsNew = copy(ds);
    dsNew.Datastore=subset(ds.Datastore,in);
    dsNew.NumObservations = numel(dsNew.Datastore.Files);
    for i = 1:dsNew.NumObservations
        file = dsNew.Datastore.Files{i};
        filepath = fileparts(file);
        [~,label] = fileparts(filepath);
        labels{i,1} = label;
    end
    dsNew.Labels = categorical(labels);
    dsNew.NumClasses = numel(unique(labels));
end

function varargout = grabEachLabel(ds,N,varargin)
    % pulls a certain number of each individual label
    opt_random=0;
    numvarargin=length(varargin);
    k=1;
    while k<=numvarargin
        switch varargin{k}
            case "random"
                opt_random=1;
        end
        k=k+1;
    end
end
```

```
    for i=1:length(N)
        in{i}=[];
    end
    labels=unique(ds.Labels);
    for i=1:ds.NumClasses
        labs=ds.Labels==labels(i);
        numLabs=sum(labs);
        x=find(labs);
        if opt_random
            x=x(randperm(length(x)));
        end
        if sum(N)>numLabs
            n=floor(numLabs.*N/sum(N));
            warning("Not enough " + string(labels(i))+" Labels")
        else
            n=N;
        end
        in{1}=[in{1}; x(1:n(1))];
        for j=2:length(n)
            in{j}=[in{j}; x(n(j-1)+1:n(j-1)+n(j))];
        end
    end
    for i=1:length(in)
        varargout{i} = subset(ds,in{i});
    end
end

function dsNew = shuffle(ds)
    % dsNew = shuffle(ds) shuffles the files and the corresponding
    % labels in the datastore.

    % Create copy of datastore.
    dsNew = copy(ds);
    dsNew.Datastore = copy(ds.Datastore);
    fds = dsNew.Datastore;

    % Shuffle files and corresponding labels.
    numObservations = dsNew.NumObservations;
    idx = randperm(numObservations);
    fds.Files = fds.Files(idx);
    dsNew.Labels = dsNew.Labels(idx);
end

end

methods (Hidden = true)
    function frac = progress(ds)
        % frac = progress(ds) returns the percentage of observations
        % read in the datastore.
```

```
        frac = (ds.CurrentFileIndex - 1) / ds.NumObservations;
    end
end

methods (Access = 'private')
    function inputStruct = parseInputs(ds,varargin)
        p = inputParser();

        p.addParameter('DataAugmentation',false,@augmentationValidator);
        p.addParameter('TestSet',false,@augmentationValidator);
        p.parse(varargin{:});
        inputStruct = p.Results;
    end
end

function data = preprocessData(Predictors,Response)
% data = preprocessData(predictors,responses) preprocesses
% the data in predictors and responses and returns the table
% data

miniBatchSize = size(Predictors,1);

% Pad data to length of longest sequence.
sequenceLengths = cellfun(@(X) size(X,2),Predictors);
maxSequenceLength = max(sequenceLengths);
for i = 1:miniBatchSize
    X = Predictors{i};

    % Pad sequence with zeros.
    if size(X,2) < maxSequenceLength
        X(:,maxSequenceLength) = 0;
    end

    Predictors{i} = X;
end

% Return data as a table.
data = table(Predictors,Response);
end

function data = readSequence(filename,newHz,aug,ts,fband,chlocs)
    % data = readSequence(filename) reads the sequence y from the MAT file
    % filename
    oldHz = 250;
    S = load(filename);
    if newHz == oldHz
        data = S.y;
    else
        data = single(resample(double(S.y),newHz,oldHz,'Dimension',2)); %✓
    end
end
```

```

end

downsample = 90; %sec
if aug
    r = randi(downsample*newHz-2)+1;
    data = data(:,r:r+(downsample*newHz-1));
end
if ts
    data = data(:,1:1+(downsample*newHz-1));
end

%Topo
data = data';
tot=bandpower(data,newHz,[1 newHz/2]);
n = size(fband,1);
abs_psd = zeros(size(data,2),n);
for i=1:n
    abs_psd(:,i)=bandpower(data,newHz,fband(i,:));
end
rel_psd = abs_psd./tot';
gscale = 67*2; %interpolation value
y = zeros(gscale,gscale,n);
for i = 1:n
    [~,y(:, :, i)] = topoplots(rel_psd(:,i),chlocs,'noplots','on','gridscale',↵
gscale);
end
data = y;
data = fillmissing(data,'constant',0);

end

function TF = augmentationValidator(valIn)

% if ischar(valIn) || isstring(valIn)
%     TF = string('none').contains(lower(valIn)); %#ok<STRQUOT>
% elseif isa(valIn,'imageDataAugmenter') && isscalar(valIn)
%     TF = true;
if islogical(valIn)
    TF = true;
else
    TF = false;
end

end
end

```